


Anti-Phospho-Rpa2 (T21) Rabbit Monoclonal Antibody For Wb, Elisa - P15927

Item Number: CM0000894



Introduction

Phospho-specific rabbit monoclonal antibody against human RPA2 phosphorylated at T21. Suitable for WB and ELISA. Cross-reacts with human and rat; detects a 29kDa band. Ideal for DNA damage and replication studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RPA2-T21 Rabbit mAb
Remarks/Alias	REPA2; RPA32; RP-A p32; RP-A p34; Phospho-RPA2-T21
Species	Human
Gene ID	6118
Immunogen	A synthetic phosphorylated peptide around T21 of human RPA2 (P15927)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P15927
Protein Weight	29kDa
Shipping	Ice bag