


Anti-Phospho-Rps6 (Ser235) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62753

Item Number: CM0001655



Introduction

Phospho-RPS6 (Ser235) rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects phosphorylated S6 in human, mouse, rat. Ideal for cell growth and proliferation studies.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-RPS6-S235 Rabbit mAb
Remarks/Alias	S6; eS6
Species	Human
GeneID (Human)	6194
GeneID	6194
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S235 of human RPS6 (NP_001001.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62753
Protein Weight	29kDa
Shipping	Ice bag