


Anti-Phospho-Rps6Ka5-S376 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O75582

Item Number: CM0001028



Introduction

Phospho-specific rabbit polyclonal antibody against RPS6KA5 (MSK1) phosphorylated at Ser376. Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying kinase signaling and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RPS6KA5-S376 Rabbit pAb
Remarks/Alias	MSK1; RLPK; MSPK1; Phospho-RPS6KA5-S376
Species	Human
GeneID (Human)	9252
GeneID	9252
Immunogen	A synthetic phosphorylated peptide around S376 of human RPS6KA5 (NP_004746.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75582
Protein Weight	90kDa
Shipping	Ice bag