

Anti-Phospho-S6 Ribosomal Protein (Rps6)-S235/236 Polyclonal Antibody For Wb, Elisa - P62753

Item Number: CM0000710



Introduction

Phospho-S6 Ribosomal Protein (RPS6)-S235/236 Rabbit pAb for WB, ELISA. Validated in human, mouse, rat. Detects phosphorylated RPS6 at Ser235/236. Ideal for cell growth and translation control studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-S6 Ribosomal Protein (RPS6)-S235/236 Rabbit pAb
Remarks/Alias	S6; eS6; Phospho-S6 Ribosomal Protein (RPS6)-S235/236
Species	Human
GeneID (Human)	6194
GeneID	6194
Immunogen	Synthetic peptide around S235 & S236 of human RPS6 (NP_001001.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62753
Protein Weight	29kDa
Shipping	Ice bag