

Anti-Phospho-S6 Ribosomal Protein (Rps6)-S235/236 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P62753

Item Number: CM0001378



Introduction

Phospho-specific rabbit monoclonal antibody recognizing human S6 Ribosomal Protein (RPS6) phosphorylated at S235/236. Validated for WB, IHC-P, and ELISA with reactivity in human, mouse, and rat. Ideal for studying cell growth and signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-S6 Ribosomal Protein (RPS6)-S235/236 Rabbit mAb
Remarks/Alias	S6, eS6, Phospho-S6 Ribosomal Protein (RPS6)-S235/236
Species	Human
GeneID (Human)	6194
GeneID	6194
Immunogen	A synthetic phosphorylated peptide around S235 & S236 of human S6 Ribosomal Protein (RPS6) (NP_001001.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62753
Protein Weight	29kDa
Shipping	Ice bag