

Anti-Phospho-S6 Ribosomal Protein (Rps6)-S240/244 Rabbit Pab For Wb, Ihc-P, Ip, Elisa - P62753

Item Number: CM0000650



Introduction

This rabbit polyclonal antibody detects endogenous levels of RPS6 only when phosphorylated at Ser240/244. Suitable for WB, IHC-P, IP, and ELISA applications in human, mouse, and rat samples. A research tool for cell growth and proliferation pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-S6 Ribosomal Protein (RPS6)-S240/244 Rabbit pAb
Remarks/Alias	S6; eS6; Phospho-S6 Ribosomal Protein (RPS6)-S240/244
Species	Human
GeneID (Human)	6194
GeneID	6194
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around S240 & S244 of human RPS6 (NP_001001.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62753
Protein Weight	29kDa
Shipping	Ice bag