


Anti-Phospho-Smad1-S463/S465+Smad5-S463/S465+Smad9-S465/S467 Rabbit Mab For Wb, Ihc-P, If/Icc, Elisa - Q15797/Q99717/O15198

Item Number: CM0001228



Introduction

Phospho-specific rabbit monoclonal antibody recognizing SMAD1 (Ser463/465), SMAD5 (Ser463/465), and SMAD9 (Ser465/467). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for BMP signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SMAD1-S463/S465+SMAD5-S463/S465+SMAD9-S465/S467 Rabbit mAb
Remarks/Alias	BSP1; MADH1; MADR1; JV4-1; DWFC; JV5-1; MADH5; MADH6; MADH9; SMAD8; Phospho-SMAD1-S463/S465+SMAD5-S463/S465+SMAD9-S465/S467
Species	Human
GenelD (Human)	4086/4090/4093
GenelD	4086/4090/4093
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15797/Q99717/O15198
Protein Weight	52kDa
Shipping	Ice bag