


Anti-Phospho-Smad2-S250 Rabbit Monoclonal Antibody For Wb, Elisa - Q15796

Item Number: CM0001721



Introduction

Phospho-Smad2-S250 Rabbit mAb validated for Western blot and ELISA. Reacts with human and mouse Smad2, phosphorylated at serine 250. Ideal for TGF-beta signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Smad2-S250 Rabbit mAb
Remarks/Alias	JV18, LDS6, CHTD8, MADH2, MADR2, JV18-1, hMAD-2, hSMAD2, Phospho-Smad2-S250
Species	Human
GeneID (Human)	4087
GeneID	4087
Immunogen	A synthetic phosphorylated peptide around S250 of human Smad2 (Q15796).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15796
Protein Weight	52kDa
Shipping	Ice bag