


Anti-Phospho-Smad2-S465/S467 Monoclonal Antibody For Wb - Q15796

Item Number: CM0001197



Introduction

Recombinant rabbit monoclonal antibody targeting phosphorylated SMAD2 at S465/S467. Validated for Western blotting (WB) and ELISA; cross-reacts with human and mouse SMAD2. Ideal for TGF-beta signaling studies.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-Smad2-S465/S467 Rabbit mAb
Remarks/Alias	JV18; LDS6; CHTD8; MADH2; MADR2; JV18-1; hMAD-2; hSMAD2; Phospho-Smad2-S465/S467
Species	Human
GeneID (Human)	4087
GeneID	4087
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S465 & S467 of human SMAD2 (NP_005892.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15796
Protein Weight	52kDa
Shipping	Ice bag