

Anti-Phospho-Smad2 (Ser467) Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15796

Item Number: CM0001062



Introduction

Rabbit polyclonal antibody against phospho-Smad2 (Ser467), validated for WB, IHC-P, and ELISA. Detects the phosphorylated TGF-beta/activin signaling transducer in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Smad2-S467 Rabbit pAb
Remarks/Alias	JV18; LDS6; CHTD8; MADH2; MADR2; JV18-1; hMAD-2; hSMAD2; Phospho-Smad2-S467
Species	Human
GeneID (Human)	4087
GeneID	4087
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S467 of human Smad2 (NP_001003652.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15796
Protein Weight	52kDa
Shipping	Ice bag