


Anti-Phospho-Smad3 (S423/S425) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P84022

Item Number: CM0000783



Introduction

Phospho-specific rabbit monoclonal antibody targeting Smad3 pS423/pS425. Validated for WB, IHC-P, and ELISA in human, mouse, and rat. Essential for studying TGF-beta/activin signaling cascades.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Smad3-S423/S425 Rabbit mAb
Remarks/Alias	LDS3; mad3; LDS1C; MADH3; JV15-2; hMAD-3; hSMAD3; HSPC193; Hst17436; Phospho-Smad3-S423/S425
Species	Human
Gene ID (Human)	4088
Gene ID	4088
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S423 & S425 of human Smad3 (P84022).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84022
Protein Weight	48 kDa
Shipping	Ice bag