


Anti-Phospho-Smad3 (S423/S425) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P84022

Item Number: CM0001513



Introduction

Phospho-specific rabbit polyclonal antibody targeting human SMAD3 phosphorylated at S423/S425. Validated for WB, IHC-P, ELISA; reacts with human, mouse, and rat. Ideal for studying TGF-beta signaling and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Smad3-S423/S425 Rabbit pAb
Remarks/Alias	LDS3; mad3; LDS1C; MADH3; JV15-2; hMAD-3; hSMAD3; HSPC193; HsT17436; Phospho-Smad3-S423/S425
Species	Human
GeneID	4088 (Human)
Immunogen	Synthetic phospho peptide corresponding to residues surrounding S423/S425 of Human SMAD3
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84022
Protein Weight	48 kDa
Shipping	Ice bag