


Anti-Phospho-Smad3-T179 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P84022

Item Number: CM0000498



Introduction

Phospho-Smad3 (Thr179) Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat Smad3. Essential for TGF-beta signaling, transcription regulation, and disease research. Reliable phospho-specific reagent.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Smad3-T179 Rabbit pAb
Remarks/Alias	LDS3; mad3; LDS1C; MADH3; JV15-2; hMAD-3; hSMAD3; HSPC193; HsT17436; Phospho-Smad3-T179
Species	Human
GeneID (Human)	4088
GeneID	4088
Immunogen	A synthetic phosphorylated peptide around T179 of human Smad3 (NP_005893.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84022
Protein Weight	48kDa
Shipping	Ice bag