


Anti-Phospho-Smad5 (S463/S465) Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99717

Item Number: CM0000787



Introduction

Rabbit polyclonal antibody specific to phosphorylated Smad5 (S463/S465). Validated for WB, IHC-P, ELISA in human, mouse, rat. Ideal for studying BMP signaling, cell differentiation, and embryonic development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Smad5-S463/S465 Rabbit pAb
Remarks/Alias	DWFC; JV5-1; MADH5; Phospho-Smad5-S463/S465
Species	Human
GeneID (Human)	4090
GeneID	4090
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S463 & S465 of human Phospho-Smad5-S463/S465 (NP_005894.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99717
Protein Weight	52kDa
Shipping	Ice bag