


Anti-Phospho-Smc1 (S957) Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q14683

Item Number: CM0001726



Introduction

Rabbit polyclonal antibody against phosphorylated SMC1A at S957, validated for WB, IHC-P, IP, ELISA. Reacts with human, mouse, rat. Suitable for cohesin, DNA repair, and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SMC1-S957 Rabbit pAb
Remarks/Alias	SMC1; SMCB; CDLS2; DEE85; SB1.8; EIEE85; SMC1L1; DXS423E; SMC1alpha; Phospho-SMC1-S957
Species	Human
Gene ID (Human)	8243
Gene ID	8243
Immunogen	A synthetic phosphorylated peptide around S957 of human SMC1A (NP_006297.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14683
Protein Weight	143kDa
Shipping	Ice bag