

Anti-Phospho-Smo (S611/S615) Rabbit Polyclonal Antibody For Wb, Elisa - Q99835

Item Number: CM0001414



Introduction

Phospho-SMO-S611/615 Rabbit Polyclonal Antibody specifically recognizes dual-phosphorylated human smoothened (SMO) at Ser611/615. Suitable for Western blot and ELISA applications in Hedgehog signaling research. Swiss-Prot: Q99835.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-SMO-S611/615 Rabbit pAb
Remarks/Alias	Gx; CRJS; PHLS; SMOH; FZD11; Phospho-SMO-S611/615
Species	Human
GeneID (Human)	6608
GeneID	6608
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99835
Protein Weight	86kDa
Shipping	Ice bag