

Anti-Phospho-Smo (S611/S615/S616) Polyclonal Antibody For Wb, Elisa - Q99835

Item Number: CM0001998



Introduction

Rabbit polyclonal antibody specific for human SMO phosphorylated at S611, S615, and S616. Validated for Western blot and ELISA, suitable for Hedgehog signaling studies. Reacts with human smoothened protein.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SMO-S611/615/616 Rabbit pAb
Remarks/Alias	Gx; CRJS; PHLS; SMOH; FZD11; Phospho-SMO-S611/615/616
Species	Human
GeneID (Human)	6608
GeneID	6608
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99835
Protein Weight	86kDa
Shipping	Ice bag