


Anti-Phospho-Sqstm1/P62-S349 Rabbit Polyclonal Antibody For Wb, Elisa - Q13501

Item Number: CM0000419



Introduction

Phospho-SQSTM1/p62-S349 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects endogenous SQSTM1/p62 phosphorylated at Ser349 in human, mouse, and rat samples. Ideal for autophagy and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SQSTM1/p62-S349 Rabbit pAb
Remarks/Alias	p60; p62; A170; DMRV; OSIL; PDB3; ZIP3; p62B; NADGP; FTDALS3; Phospho-SQSTM1/p62-S349
Species	Human
Gene ID (Human)	8878
Immunogen	A synthetic phosphorylated peptide around S349 of human SQSTM1/p62 (NP_003891.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13501
Protein Weight	48kDa
Shipping	Ice bag