


Anti-Phospho-Src Family (Y416) Rabbit Monoclonal Antibody For Wb, Elisa - P12931 / P06241 / P06239 / P07947 / P08631 / P07948

Item Number: CM0000811



Introduction

Phospho-specific rabbit monoclonal antibody detecting SRC family kinases at the conserved Y416 activation loop (human, cross-reacts with rat). Validated for western blot and ELISA, supporting signal transduction and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SRC Family-Y416 Rabbit mAb
Remarks/Alias	ASV; SRC1; THC6; c-SRC; p60-Src; SRC; Phospho-SRC Family-Y416
Species	Human
Gene ID (Human)	2534/3055/3932/4067/6714/7525
Gene ID	2534/3055/3932/4067/6714/7525
Immunogen	A synthetic phosphorylated peptide around Y416 of human SRC Family (NP_005408.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P12931/P06241/P06239/P07947/P08631/P07948
Protein Weight	60kDa
Shipping	Ice bag