

Anti-Phospho-Src (S17) Rabbit Polyclonal Antibody For Wb And Elisa - P12931

Item Number: CM0001886



Introduction

Phospho-Src-S17 Rabbit Polyclonal Antibody validated for WB and ELISA in human, mouse, rat. Ideal for studying Src kinase signaling in cell adhesion, cancer, and immune response.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Src-S17 Rabbit pAb
Remarks/Alias	ASV; SRC1; THC6; c-SRC; p60-Src; Phospho-Src-S17
Species	Human
GeneID (human)	6714
GeneID	6714
Immunogen	A synthetic phosphorylated peptide around S17 of human Src (NP_005408.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12931
Protein Weight	60kDa
Shipping	Ice bag