

Anti-Phospho-Src-Y529 Rabbit Monoclonal Antibody For Wb And Elisa - P12931

Item Number: CM0000793



Introduction

High-quality rabbit monoclonal antibody against phosphorylated Src at tyrosine 529 for Western blot and ELISA. Reacts with human samples. Ideal for cancer, cell adhesion, and bone resorption research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Src-Y529 Rabbit mAb
Remarks/Alias	ASV; SRC1; THC6; c-SRC; p60-Src; Phospho-Src-Y529
Species	Human
GeneID (Human)	6714
GeneID	6714
Immunogen	Synthetic Peptide - A synthetic phosphorylated peptide around Y529 of human Src.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P12931
Protein Weight	60kDa
Shipping	Ice bag