

Anti-Phospho-Stat1 (S727) Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P42224

Item Number: CM0001317



Introduction

Phospho-specific rabbit monoclonal antibody targeting STAT1 phosphorylated at Ser727. Validated for WB, IF/ICC, IHC-P, ELISA in human, mouse, rat. Ideal for IFN signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT1-S727 Rabbit mAb
Remarks/Alias	CANDF7; IMD31A; IMD31B; IMD31C; ISGF-3; STAT91; Phospho-STAT1-S727
Species	Human
GeneID (human)	6772
GeneID	6772
Immunogen	A synthetic phosphorylated peptide around S727 of human STAT1 (P42224)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42224
Protein Weight	87kDa
Shipping	Ice bag