


Anti-Phospho-Stat1 (Ser727) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P42224

Item Number: CM0001757



Introduction

Rabbit polyclonal antibody targeting phosphorylated STAT1 (Ser727) for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat samples. Expected molecular weight 87 kDa. UniProt P42224.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT1-S727 Rabbit pAb
Remarks/Alias	CANDF7; IMD31A; IMD31B; IMD31C; ISGF-3; STAT91; Phospho-STAT1-S727
Species	Human
Gene ID (Human)	6772
Gene ID	6772
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S727 of human STAT1 (NP_009330.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42224
Protein Weight	87 kDa
Shipping	Ice bag