

Anti-Phospho-Stat1 (Y701) Rabbit Polyclonal Antibody For Wb, IHC-P, Elisa - P42224

Item Number: CM0000961



Introduction

Phospho-STAT1 (Y701) rabbit polyclonal antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Detects phosphorylated STAT1, a key transcription factor in interferon signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT1-Y701 Rabbit pAb
Remarks/Alias	CANDF7; IMD31A; IMD31B; IMD31C; ISGF-3; STAT91; Phospho-STAT1-Y701
Species	Human
GeneID (human)	6772
GeneID	6772
Immunogen	A synthetic phosphorylated peptide around Y701 of human STAT1 (NP_009330.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42224
Protein Weight	87kDa
Shipping	Ice bag