

Anti-Phospho-Stat2 (Tyr690) Rabbit Monoclonal Antibody For Wb And Elisa - P52630

Item Number: CM0001620



Introduction

Phospho-specific rabbit monoclonal antibody against STAT2 (Tyr690) for WB and ELISA. Validated in human and rat samples; detects 98 kDa band. Ideal for type I interferon signaling and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT2-Y690 Rabbit mAb
Remarks/Alias	P113, IMD44, ISGF-3, PTORCH3, STAT113
Species	Human
GeneID (Human)	6773
GeneID	6773
Immunogen	A synthetic phosphorylated peptide around Y690 of human STAT2 (NP_005410.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P52630
Protein Weight	98 kDa
Shipping	Ice bag