

Anti-Phospho-Stat3-S727 Rabbit Mab - P40763

Item Number: CM0000868



Introduction

Phosphorylated STAT3 (S727) rabbit monoclonal antibody validated for WB, IP, IF/ICC, IHC-P, ChIP, ELISA. Detects human, mouse, rat STAT3. Ideal for JAK/STAT signaling, inflammation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT3-S727 Rabbit mAb
Remarks/Alias	APRF; HIES; ADMIO; ADMIO1
Species	Human
GeneID (human)	6774
GeneID	6774
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S727 of human STAT3 (NP_644805.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, IF/ICC, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40763
Protein Weight	88 kDa/83 kDa
Shipping	Ice bag