

Anti-Phospho-Stat3-S727 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P40763

Item Number: CM0001553



Introduction

Rabbit monoclonal antibody specific for STAT3 phosphorylated at Ser727. Validated for WB, IHC-P, ELISA with human, mouse, rat reactivity. Ideal for investigating STAT3 signaling in inflammation, cell cycle, and T-cell differentiation.

[Learn More](#)

Parameter	Description
Product Name	Phospho-STAT3-S727 Rabbit mAb
Remarks/Alias	APRF; Phospho-STAT3-S727
Species	Human
Gene ID	6774
Immunogen	A synthetic phosphorylated peptide around S727 of human STAT3 (P40763)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40763
Protein Weight	88 kDa
Shipping	Ice bag