


Anti-Phospho-Stat5-Y694 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P42229/P51692

Item Number: CM0001405



Introduction

Phospho-STAT5-Y694 Rabbit mAb detects human STAT5A/B phosphorylated at Y694. Validated for WB, IF-P, IHC-P, ELISA. For studying JAK/STAT pathway activation and transcriptional regulation.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-STAT5-Y694 Rabbit mAb
Remarks/Alias	MGF; STAT5; stat5a/stat5b; Phospho-STAT5-Y694
Species	Human
GeneID (human)	6776/6777
GeneID	6776/6777
Immunogen	A synthetic phosphorylated peptide around Y694 of human STAT5 (NP_001275647.1/NP_036580.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P42229/P51692
Protein Weight	90kDa
Shipping	Ice bag