


Anti-Phospho-Stat5-Y694 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P42229/P51692

Item Number: CM0001026



Introduction

Phospho-STAT5-Y694 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Detects phosphorylated STAT5A/STAT5B at Y694, key for cytokine signaling and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT5-Y694 Rabbit pAb
Remarks/Alias	MGF; STAT5; stat5a/stat5b; Phospho-STAT5-Y694
Species	Human
GeneID (Human)	6776/6777
GeneID	6776/6777
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y694 of human STAT5 (NP_003143.2/NP_036580.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42229/P51692
Protein Weight	90kDa
Shipping	Ice bag