


Anti-Phospho-Stat5A (S780) Rabbit Polyclonal Antibody For Wb - P42229

Item Number: CM0000630



Introduction

Phospho-STAT5A (S780) rabbit polyclonal antibody validated for Western blot in human samples. Detects the activated form of STAT5A, a key transcription factor in cytokine signaling and lactation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT5A-S780 Rabbit pAb
Remarks/Alias	MGF; STAT5; Phospho-STAT5A-S780
Species	Human
Gene ID (Human)	6776
Gene ID	6776
Immunogen	Synthetic phospho-specific peptide corresponding to residues surrounding S780 of human STAT5A
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P42229
Protein Weight	91kDa
Shipping	Ice bag