

Anti-Phospho-Stathmin 1 (S25) Rabbit Monoclonal Antibody For Wb, Elisa - P16949

Item Number: CM0001203

Introduction

Rabbit monoclonal phospho-specific antibody targeting Stathmin 1 (S25) for WB and ELISA, recognizing human and rat. Suitable for studying microtubule dynamics, cell cycle regulation, and neurogenesis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Stathmin 1-S25 Rabbit mAb
Remarks/Alias	Lag; SMN; OP18; PP17; PP19; PR22; LAP18; C1orf215; Phospho-Stathmin 1-S25
Species	Human
GeneID (Human)	3925
GeneID	3925
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S25 of human Stathmin 1 (NP_005554.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P16949
Protein Weight	17kDa
Shipping	Ice bag