

Anti-Phospho-Stathmin 1 (Ser16) Rabbit Polyclonal Antibody For Wb And Elisa - P16949

Item Number: CM0000620

Introduction

Phospho-Stathmin 1 (Ser16) rabbit polyclonal antibody validated for western blot (WB) and ELISA. Detects human and mouse Stathmin 1, a microtubule-destabilizing phosphoprotein critical for neurogenesis and frequently overexpressed in leukemia.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Stathmin 1-S16 Rabbit pAb
Remarks/Alias	Lag; SMN; OP18; PP17; PP19; PR22; LAP18; C1orf215; Phospho-Stathmin 1-S16
Species	Human
GeneID (Human)	3925
GeneID	3925
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S16 of human Stathmin 1 (NP_001138926.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16949
Protein Weight	17kDa
Shipping	Ice bag