

Anti-Phospho-Stathmin 1 (Ser25) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P16949

Item Number: CM0001090

Introduction

Validated for Western blot, IHC-P, and ELISA, this rabbit polyclonal antibody detects endogenous phospho-Stathmin 1 (Ser25) in human, mouse, and rat samples, supporting cancer and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Stathmin 1-S25 Rabbit pAb
Remarks/Alias	Lag; SMN; OP18; PP17; PP19; PR22; LAP18; C1orf215; Phospho-Stathmin 1-S25
Species	Human
Gene ID (Human)	3925
Gene ID	3925
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S25 of human Stathmin 1 (NP_005554.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16949
Protein Weight	17kDa
Shipping	Ice bag