

Anti-Phospho-Sting1 (Ps366) Rabbit Monoclonal Antibody For Wb, Elisa - Q86Wv6

Item Number: CM0001989



Introduction

Phospho-STING/TMEM173-S366 Rabbit mAb detects endogenous STING1 protein when phosphorylated at Ser366. Suitable for WB and ELISA, this antibody is a critical tool for studying innate immune signaling, cGAS-STING pathway, and host-pathogen interactions.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STING/TMEM173-S366 Rabbit mAb
Remarks/Alias	ERIS; MITA; MPYS; SAVI; NET23; STING; hMITA; hSTING; TMEM173; STING-beta; Phospho-STING/TMEM173-S366
Species	Human
Gene ID (Human)	340061
Gene ID	340061
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S366 of human STING1 (NP_938023.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86WV6
Protein Weight	42kDa
Shipping	Ice bag