


Anti-Phospho-Stk11-T189 Rabbit Pab For Wb - Q15831

Item Number: CM0001659



Introduction

Rabbit polyclonal antibody detecting phosphorylation at T189 of human STK11 (LKB1), validated for WB and ELISA. Reacts with human, mouse, rat. Ideal for AMPK signaling, tumor suppression, and cell polarity studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STK11-T189 Rabbit pAb
Remarks/Alias	PJS; LKB1; hLKB1; Phospho-STK11-T189
Species	Human
Gene ID (Human)	6794
Gene ID	6794
Immunogen	A synthetic phosphorylated peptide around T189 of human STK11 (NP_000446.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15831
Protein Weight	49kDa
Shipping	Ice bag