

Anti-Phospho-Stk3/Stk4-T180/T183 Polyclonal Antibody For Wb And Elisa - Q13188 / Q13043

Item Number: CM0001696



Introduction

High-affinity rabbit polyclonal antibody specific for phosphorylated STK3/STK4 at T180/T183, validated for WB and ELISA in human, mouse, and rat samples. Ideal for studying Hippo signaling, apoptosis, and centrosome biology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STK3/STK4-T180/T183 Rabbit pAb
Remarks/Alias	KRS1; MST2; KRS2; MST1; YSK3
Species	Human
GeneID (Human)	6788/6789
GeneID	6788/6789
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T180 of human STK3/STK4STK3 (NP_006272.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13188/Q13043
Protein Weight	56kDa
Shipping	Ice bag