

Anti-Phospho-Syk-Y323 Polyclonal Antibody For Wb, Elisa - P43405

Item Number: CM0001609



Introduction

Rabbit polyclonal antibody specific to human Syk phosphorylated at Tyr323. Validated for Western blot and ELISA; detects a 72 kDa band. Does not cross-react with other species. Suitable for B-cell signaling and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Syk-Y323 Rabbit pAb
Remarks/Alias	IMD82; p72-Syk; Phospho-Syk-Y323
Species	Human
GeneID (Human)	6850
GeneID	6850
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y323 of human Syk (NP_003168.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P43405
Protein Weight	72kDa
Shipping	Ice bag