

Anti-Phospho-Syk (Y352)/Zap70 (Y319) Rabbit Polyclonal Antibody (Pab) For Wb, Elisa - P43405 / P43403

Item Number: CM0001843



Introduction

Rabbit polyclonal antibody targeting phosphorylated Syk (Y352) and ZAP70 (Y319) for Western blot and ELISA applications. Specific for human samples, ideal for immunoreceptor signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Syk-Y352/ZAP70-Y319 Rabbit pAb
Remarks/Alias	SYK/ZAP70; Phospho-Syk-Y352/ZAP70-Y319
Species	Human
Gene ID (Human)	6850/7535
Gene ID	6850/7535
Immunogen	A synthetic phosphorylated peptide around Y352 of human Syk (NP_003168.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P43405/P43403
Protein Weight	69kDa/72kDa/35kDa/55kDa
Shipping	Ice bag