

Anti-Phospho-Syk (Y525/526) Rabbit Polyclonal Antibody For Wb And Elisa - P43405

Item Number: CM0000838



Introduction

Phospho-specific rabbit polyclonal antibody against human SYK phosphorylated at Tyr525/526. Validated for Western blot and ELISA. Ideal for B-cell receptor signaling, innate immunity, and kinase activity studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Syk-Y525/526 Rabbit pAb
Remarks/Alias	IMD82; p72-Syk; Phospho-Syk-Y525/526
Species	Human
GeneID (human)	6850
GeneID	6850
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y525 & Y526 of human Syk (NP_003168.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P43405
Protein Weight	72kDa
Shipping	Ice bag