

Anti-Phospho-Synapsin-1 (Ser9) Polyclonal Antibody For Wb And Elisa - P17600

Item Number: CM0001262

Introduction

High-quality rabbit polyclonal antibody targeting phosphorylated Synapsin-1 at Ser9. Validated for WB and ELISA, recognizing human, mouse, and rat species. Ideal for synaptic vesicle trafficking and neurotransmitter release studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SYN1-S9 Rabbit pAb
Remarks/Alias	SYNI; EPILX; MRX50; SYN1a; SYN1b; EPILX1; Phospho-SYN1-S9
Species	Human
GeneID (human)	6853
GeneID	6853
Immunogen	Synthetic peptide—A synthetic phosphorylated peptide around S9 of human SYN1 (NP_008881.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17600
Protein Weight	74kDa
Shipping	Ice bag