


Anti-Phospho-Tak1-S439 Rabbit Monoclonal Antibody For Wb - O43318

Item Number: CM0000553



Introduction

Phospho-TAK1-S439 rabbit monoclonal antibody validated for Western blot and ELISA, cross-reacts with mouse. Designed for studying MAPK signaling, NF- κ B activation, and cellular stress responses.

[Learn More](#)

Core Parameter	Value
Product Name	Phospho-TAK1-S439 Rabbit mAb
Remarks/Alias	CSCF; FMD2; TAK1; MEKK7; TGF1a; Phospho-TAK1-S439
Species	Human
GeneID (Human)	6885
GeneID	6885
Immunogen	A synthetic phosphorylated peptide around S439 of human TAK1 (O43318)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O43318
Protein Weight	67kDa
Shipping	Ice bag