

Anti-Phospho-Tau (S396) Rabbit Monoclonal Antibody For Wb, Elisa - P10636

Item Number: CM0001451



Introduction

Phospho-Tau (S396) rabbit monoclonal antibody for detection of phosphorylated tau protein, validated in WB and ELISA. Cross-reacts with mouse, suitable for neurodegeneration research and assay development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-S396 Rabbit mAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-S396
Species	Human
GeneID (Human)	4137
GeneID	4137
Immunogen	A synthetic phosphorylated peptide around S396 of human tau (NP_005901.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag