

Anti-Phospho-Tau-S404 Rabbit Monoclonal Antibody For Wb And Elisa - P10636

Item Number: CM0001446



Introduction

High-quality Phospho-Tau-S404 Rabbit mAb validated for WB and ELISA. Recognizes human, mouse, rat phosphorylated Tau at S404. Ideal for neurodegenerative disease research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-S404 Rabbit mAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-S404
Species	Human
Gene ID (Human)	4137
Immunogen	A synthetic phosphorylated peptide around S404 of human Tau
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag