

Anti-Phospho-Tau (Ser416) Rabbit Polyclonal Antibody For Wb, Ihc, If, Elisa - P10636

Item Number: CM0001732



Introduction

Rabbit polyclonal antibody specific for human Tau phosphorylated at Ser416. Validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Detects human, mouse, and rat MAPT. For neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-Ser416 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-Ser416
Species	Human
GeneID (Human)	4137
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S416 of human Tau (NP_005901.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag