

Anti-Phospho-Tau (T205) Monoclonal Antibody For Wb, Ihc-P, If - P10636

Item Number: CM0001241



Introduction

Rabbit monoclonal antibody against phosphorylated Tau (T205). Validated for WB, IHC-P, IF, and ELISA. Reacts with human, mouse, and rat. Ideal for Alzheimer's disease and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T205 Rabbit mAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T205
Species	Human
GeneID (Human)	4137
GeneID	4137
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T205 of human Tau (NP_001190181.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P10636
Protein Weight	79 kDa
Shipping	Ice bag