

Anti-Phospho-Tau (T205) Rabbit Polyclonal Antibody For Wb, Elisa - P10636

Item Number: CM0001678



Introduction

Phospho-Tau (T205) rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Detects human, mouse, and rat tau phosphorylated at threonine 205. Ideal for Alzheimer's and tauopathy research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T205 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T205
Species	Human
GeneID (human)	4137
GeneID	4137
Immunogen	A synthetic phosphorylated peptide around T205 of human Tau (NP_001190181.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag