

Anti-Phospho-Tau-T212 Rabbit Polyclonal Antibody For Wb And Elisa - P10636

Item Number: CM0001718



Introduction

Rabbit polyclonal antibody targeting human Tau phosphorylated at T212. Validated for WB and ELISA with cross-reactivity to mouse and rat. Useful for studying microtubule dynamics, neuronal polarity, and tauopathies including Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T212 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T212
Species	Human
GeneID	4137
Immunogen	Synthetic peptide - A phospho specific peptide corresponding to residues surrounding T212 of human Tau
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag