

Anti-Phospho-Tau-T217 Rabbit Polyclonal Antibody (Pab) For Wb, Elisa - P10636-5

Item Number: CM0001504

Introduction

Phospho-Tau-T217 Rabbit Polyclonal Antibody (pAb) for WB, ELISA. Detects phosphorylated Tau at T217 in human, mouse, and rat samples. Ideal for Alzheimer's and tauopathy research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T217 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T217
Species	Human
GeneID (Human)	4137
GeneID	4137
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around T217 of human Tau (NP_001190181.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636-5
Protein Weight	32-45kDa/78-80kDa
Shipping	Ice bag