

Anti-Phospho-Tau-T231 Monoclonal Antibody For If-P, Elisa - P10636

Item Number: CM0000788



Introduction

Phospho-Tau-T231 Rabbit monoclonal antibody validated for IF-P and ELISA with human, mouse, rat reactivity. Targets phosphorylated Thr231 of Tau (P10636), key in Alzheimer's and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T231 Rabbit mAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T231
Species	Human
Gene ID	4137
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around T231 of human Tau (P10636).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag