

Anti-Phospho-Tau-Thr212/Ser214 Rabbit Polyclonal Antibody For Wb, Elisa - P10636

Item Number: CM0000467



Introduction

Rabbit polyclonal antibody against phosphorylated human Tau (Thr212/Ser214), validated for WB and ELISA. Cross-reacts with mouse Tau. Ideal for Alzheimer's and neurodegeneration research. Protein weight ~79 kDa.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-Thr212/Ser214 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-Thr212/Ser214
Species	Human
GeneID	4137 (Human)
Immunogen	Synthetic phosphorylated peptide around T529 & S531 of human Tau (NP_058519.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10636
Protein Weight	79 kDa
Shipping	Ice bag